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



STERRAD® Plasma Sterilization

ASP – the inventor of low temperature hydrogen peroxide plasma sterilization technology

STERRAD Sterilization – worldwide safe, proven and experienced plasma technology

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SPECTRUM OF ACTIVITY

TABLE 1
VEGETATIVE BACTERIAL, SPORES AND FUNGI

MICROORGANISM	TYPE	INTEREST IN TESTING	CONTROL*	RESULTS**
<i>Geobacillus stearothermophilus</i>	Bacterial Spore	H ₂ O ₂ Resistance; Stream Indicator Organism	2.04 x 10 ⁶	0/9
<i>Bacillus subtilis</i> var. <i>niger</i> (globigii)	Bacterial Spore	H ₂ O ₂ Resistance; EO Indicator Organism	2.69 x 10 ⁶	0/9
<i>Bacillus pumilus</i>	Bacterial Spore	Ionizing Radiation Resistance; Indicator Organism	1.82 x 10 ⁶	0/9
<i>Staphylococcus aureus</i>	Gram Positive	H ₂ O ₂ Resistance; Clinical Significance	2.82 x 10 ⁶	0/9
<i>Deinococcus radiodurans</i>	Gram Positive	Ionizing Radiation Resistance	3.10 x 10 ⁶	0/9
<i>Pseudomonas aeruginosa</i>	Gram Negative	Clinical Significance	1.32 x 10 ⁶	0/9
<i>Escherichia coli</i>	Gram Negative	Clinical Significance	9.23 x 10 ⁶	0/9
<i>Serratia marcescens</i>	Gram Negative	H ₂ O ₂ Resistance; Clinical Significance	1.85 x 10 ⁶	0/9
<i>Moraxella osloensis</i>	Gram Negative	Ionizing Radiation Resistance EO Indicator Organism	3.14 x 10 ⁶	0/9
<i>Mycobacterium bovis</i>	Acid Fast	Chemical Resistance; Clinical Resistance	4.20 x 10 ⁶	0/9
<i>Candida albicans</i>	Yeast	H ₂ O ₂ Resistance	3.95 x 10 ⁶	0/9
<i>Candida parapsilosis</i>	Yeast	H ₂ O ₂ Resistance; Clinical Significance	1.07 x 10 ⁶	0/9
<i>Trichophyton mentagrophytes</i>	Filamentous Fungus	Clinical Significance	1.25 x 10 ⁶	0/9
<i>Aspergillus niger</i>	Filamentous Fungus	H ₂ O ₂ Resistance; Clinical Significance	1.46 x 10 ⁶	0/9

* Average titer recovered from nine samples.
** # Positive/ # Tested.

TABLE 1 (CONTINUED)
VIRUSES

MICROORGANISM	TYPE	INTEREST IN TESTING	VIRUS TITER LOG 10	INFECTIVITY
<i>Poliovirus Type 1</i> (Brunhilde)	Hydrophilic	Chemical Resistance; Clinical Significance	Test 1 >3.98 Test 2 >3.98	Not detected Not detected
<i>Herpesvirus Type 1</i>	Lipophilic	Clinical Significance	Test 1 >3.20 Test 2 >2.84	Not detected Not detected

TABLE 2
BACTERIA

Clinical Isolates Tested in an Abbreviated STERRAD® System Cycle

MICROORGANISM	#POSITIVE/#TESTED**	MICROORGANISM	#POSITIVE/#TESTED**
<i>Pseudomonas aeruginosa</i>	0/10	<i>Aeromonas</i> sp.	0/10
<i>Pseudomonas cepacia</i>	0/10	<i>Clostridium perfringens</i>	0/10
<i>Xanthomonas maltophilia</i>	0/10	<i>Clostridium tetani</i>	0/10
<i>Serratia marcescens</i>	0/10	<i>Clostridium difficile</i>	0/10
<i>Klebsiella (encapsulated)</i>	0/10	<i>Bacillus subtilis</i> spores	0/10
<i>Methicillin-resistant Staphylococcus aureus</i>	0/10	<i>Micrococcus</i> sp.	0/10
<i>Slime-producing Staphylococcus epidermidis</i>	0/10	<i>Mycobacterium tuberculosis</i>	0/10
<i>Listeria monocytogenes</i>	0/10	<i>Mycobacterium chelonae</i>	0/10
<i>Enterococcus faecalis</i>	0/10	<i>Bacteroides fragilis</i>	0/10
<i>Acinetobacter calcoaceticus</i>	0/10	<i>Fusobacterium</i> sp.	0/10
<i>Salmonella</i> sp.	0/10	<i>Anaerobic cocci</i>	0/10
<i>Shigella</i> sp.	0/10	<i>Candida albicans</i>	0/10
<i>Campylobacter</i> sp.	0/10		

* Abbreviated cycle consisted of 40 minutes of diffusion with 3 mg/liter of hydrogen peroxide and 10 minutes of plasma at a power of 300 watts.
** Five tests were conducted in the presence of 5% serum and five tests were conducted in the presence of 10% serum.

TABLE 3 A

Half-Cycle Validation Test with 1.3×10^6 *G. stearothermophilus* Spores on Stainless Steel Carriers in 0.7 mm x 500 mm Stainless Steel Lumens in a Standardized Validation Load with Reduced Concentrations of Hydrogen Peroxide.

STERRAD® 100NX™ SYSTEM STANDARD CYCLE

Half-Cycle Number	ml of H2O2 Injection 53%	Enumeration Test Spore Log Reduction	Biological Results # non-sterile / # tested
1	06.2	1.39	-
2	0.93	3.94	-
3	1.23	-	5/10
4	2.45	-	0/10
5	3.10	-	0/10
6	3.68	-	0/10
7	4.26	-	0/10
8	4.90	-	0/10

TABLE 3 B

Half-Cycle Validation Test with 2.6×10^6 *G. stearothermophilus* Spores on Stainless Steel Carriers in 1 mm x 1000 mm Polyethylene and Teflon Tubing with Reduced Concentrations of Hydrogen Peroxide.

STERRAD® 100NX™ SYSTEM STANDARD CYCLE

Half-Cycle Number	ml of H ₂ O ₂ Injection 53%	Biological Results # non-sterile / # tested
1	0.93	0/10
2	1.23	0/10
3	2.45	0/10
4	3.10	0/10
5	3.68	0/10
6	4.26	0/10
7	4.90	0/20

TABLE 4

G. stearothermophilus Spores between Mated Surfaces in a Standardized Validation Load at Half-Cycle Conditions

Material	Test		
	# Positive / # Tested	# Positive / # Tested	# Positive / # Tested
Stainless Steel Blades	0/10	0/10	0/10
Titanium	0/10	0/10	0/10
Polyacetal (Delrin)	0/10	0/10	0/10
Polyetherimide (Ultem*)	0/10	0/10	0/10
Polyarylsulfone (Radel**)	0/10	0/10	0/10

*Ultem is a trademark of GE Plastics **Radel is a trademark of Amoco

TABLE 5

Half-Cycle Validation Test with 2.2×10^6 *G. stearothermophilus* Spores on Stainless Steel Carriers in 1 mm x 850 mm Flexible Endoscope with Reduced Concentrations of Hydrogen Peroxide.

STERRAD® 100NX™ SYSTEM FLEX SCOPE CYCLE

Cycle # (per device)	Biological Results: # non-sterile / # Tested			
	Choledochoscope	Ureteroscope	Bronchovideoscope	Ureterorenofiberscope
1	0/1	0/1	0/2	0/1
2	0/1	0/1	0/2	0/1
3	0/1	0/1	0/2	0/1

TABLE 6

AOAC Sporicidal Test with STERRAD® 100NX™ System

Test #	Carrier	# of Failures / # Tested	
		<i>B. Subtilis</i>	<i>C. Sporogenes</i>
1	Suture	0/60	0/60
	Porcelain	0/60	0/60
2	Suture	0/60	0/60
	Porcelain	0/60	0/60
3	Suture	0/60	0/60
	Porcelain	0/60	0/60

Medical Device Sterility Assurance

- ASP's Medical Device Manufacturer Team - mdmprogram@its.jnj.com
- MDM Testing Program - mdmprogram@its.jnj.com
- STERRAD® Sterility Guide (SSG)- www.sterradsterilityguide.com

1. Instrument preparation

- **Disassemble** instruments according to the manufacturer's Instructions for Use (IFU). **Remove** all blood, tissue and soil from instruments using an **enzymatic detergent** per the manufacturer's IFU. **Rinse** the instruments thoroughly to remove any residue.
- **Dry all items completely** and inspect each instrument for proper function. Moisture on instruments will lead to cancellation of the cycle. **Medical grade compressed air** may be used to help dry lumens and other hidden spaces.


2. Instrument packaging: I


- Instruments can be **packaged** for sterilization utilizing different methods.
- **APTIMAX® Instrument Trays** are recommended for use in the STERRAD® System, as they have been designed to allow optimal diffusion of hydrogen peroxide gas around all items in the load. Use APTIMAX® Instrument Tray Holders to secure instruments in the trays. APTIMAX® Instrument Trays can be enclosed in either **polypropylene wrap** or in **Tyvek® Pouches** with STERRAD® **Chemical Indicators**.



Arthroscopy set

Cystoscopy set

Hysteroscopy set

2. Instrument packaging: II

- Container systems tested and cleared for use in STERRAD® Systems, such as ASP SteriTite® Containers,* may be useful for processing typical instrument sets.
- Place a STERRAD® Chemical Indicator (CI) Strip inside each container system or tray prior to wrapping to provide confirmation that the instrument has been processed.

2. Instrument packaging: III

- If the instrument tray has been wrapped in polypropylene wrap, **secure the wrap with STERRAD® SEALSURE® Tape**. The tape color will **change from red to gold** Chemical Indicator (CI) or lighter during the sterilization cycle, indicating that the tray has been exposed to hydrogen peroxide.
- **Single instruments may be placed in peel pouches**, enabling easy identification of pouch contents. **Peel pouches should not be placed inside wrapped trays**.

2. Instrument packaging: IV

- **Tyvek® Pouches** with STERRAD® Chemical Indicators are made of Tyvek®, which **is permeable to hydrogen peroxide vapor** and provides a microbial barrier.
- A STERRAD® Chemical Indicator Strip should be placed **inside each pouch**.
- **Writing the contents** on the STERRAD® SEALSURE®
- Chemical Indicator (CI) Tape may assist with identification of the pouch contents.

3. Loading the sterilizer

- Multiple **peel pouches** should be placed **on edge inside an uncovered tray**. Pouches should be placed in the tray so that the clear side of one pouch faces the opaque side of the next pouch. A single-peel pouch should be placed flat on the shelf with the Tyvek® side facing up.
- Provide at least **1 inch (2.5 mm) of space** between the electrode and the top of the load, 1 inch between packages in the load and 1 inch between the electrode and the sides of the load. Place trays flat on the sterilizer shelves in a single layer. Do not stack trays.

Do not process the following

- Items made of materials containing **cellulose**, such as cotton, paper or cardboard, linens, huck towels, gauze sponges, or any item containing wood pulp
- **Single-use items**
- Instruments with lumens that have **not been completely dried**
- Items wrapped in barrier material containing **paper**
- **Foam pads** in instrument trays

When loading the sterilizer, DO NOT :

- **Stack** instrument trays
- Place instrument trays **on edge** in order to fit more items in the sterilizer
- Insert the sterilizer shelf **upside-down**
- Configure the load so that metal instruments are **in contact with the walls or door of the sterilizer**

TABLE 7


Representative Medical Devices Evaluated for Compatibility with the STERRAD® 100NX™ Sterilization System after Repetitive Sterilizations

DEVICE	PROPERTY MEASURED	TOTAL CYCLES
Resectoscope	Power Output	50
Forceps	Mechanical Properties	50
Defibrillator Handle and Internal Electrodes	Power Output Compared to Charge Taken	50
Microsurgical Instruments	Subjective Evaluation of Instrument Appearance and Sharpness by Visual Observation, and Mechanical Function of Moving Parts	50
Fiberoptic Ureteroscope	Optics, Mechanical Properties, Subjective Evaluation of Instruments Appearance	50
Rigid Hysteroscope	Optics, Mechanical Properties, Subjective Evaluation of Instruments Appearance	50

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Sterility Assurance Levels


FIGURE 1

Graphic Presentation of Half-Cycle Method of Sterilizer (2,3)



FIGURE 2

Graphic Presentation of the Two-Phase Cycle

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
對殺菌化學劑抵抗性的遞減順序 (Descending Order of Resistance to Germicidal Chemicals)

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Network Connectivity

- Customers have the option to connect their STERRAD 100NX System to a remote PC or network via an Ethernet connection on the sterilizer.
- Ability to electronically maintain cycle records and download & view one-second cycle data.
- Allows for easy upload of instrument inventory to create menu for load entry data list. Eliminates need to hand-write load list.

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- An Optional IMS is available.
 - An IMS is a combination of process measurement instrumentation, software and recording device that duplicates the sterilizer control system instrumentation.
 - 獨立監測系統結合了程序測量設備、軟體及記錄裝置，以便複製滅菌器的控制系統設備資訊。
 - There are 9 sensors included in the IMS kit.
 - 3 sensors monitor pressure, 5 monitor temperature and 1 monitors plasma power.
 - IMS 套件中共有 9 個感應器，其中 3 個感應器監測壓力，5 個感應器監測溫度，並有 1 個感應器監測電漿電源。

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Typical Devices Sterilized in the STERRAD® 100NX™ Sterilizer*

- Cranial pressure transducer cables
- Cryoprobes
- Defibrillator paddles
- Dopplers
- Electrocautery instruments
- Endoscopic instruments
- Esophageal dilators
- Fiberoptic light cables
- Laryngoscope blades
- Laser hand-pieces, fibers, and accessories
- Metal instruments
- Ophthalmic lenses (diagnostic, magnifying)
- Patient lead cables
- Pigmentation hand-pieces
- Radiation therapy equipment
- Resectoscope/working elements and sheaths
- Rigid endoscopes
- Shaver hand-pieces
- Single-channel flexible endoscopes
- Stereotactic equipment
- Surgical power equipment and batteries
- Trocar sheaths
- Ultrasound probes
- Video cameras and couplers

*Any devices processed in the STERRAD® 100NX™ Sterilizer must be within the claim limits of the sterilizer.

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